

eatures – Web Site User Account Login et Fall Detection & Fall Detection Sensitivity	Basic	Plus
at Fall Detection & Fall Detection Constituity		rius
et Fail Detection & Fail Detection Sensitivity	\checkmark	✓
et location upload interval (default set at 60 minutes)	✓	✓
et Geo-Fence and receive alerts if Wristband Wearer wanders out of Geo-Fence area.	✓	✓
et up to 3 Medication Reminders by time of day & day of week.	✓	✓
dd Family & Friends phone numbers & email. Receive confirmation of acceptance of neir participation in your Emergency Contacts group.	✓	✓
Allow a limited number of Family & Friends to call the Wristband directly.		✓
end email invitation to Family & Friends to download the App.	√	✓
iew and print out a personalized user guide detailing how the Wristband is set up with dvice to get the most out of the Wristband.	✓	✓
iew and print out Credit-Card sized User Guide for Friends & Family members.		✓
iew connectivity statistics for Wristband mobile signal level.		✓
iew battery charge statistics with highlights of any charging problems.		✓
iew activity level statistics.		✓
	Flawless	Flawless
mergency Call Dial Out Features	Basic	Plus
apid dial out to multiple callers in turn (approx. 45 seconds to each contact in turn).	√ V	√
mergency Contact presses 1 to avoid call going to Voice Mail.	√	✓
all All & Connect To First – fastest contact method.		√
SMS sent to Family & Friends who received call but didn't answer confirming Wristband		✓
was connected to Emergency Contact OK including which contact answered. Rapid bring in Emergency Services (999,000,911) on 3-way call.		✓
Rapid bring in other Emergency Contacts on 3-way call via short code dial.		✓
apid precise instant location detection using What3Words using advanced text-to- peech over the phone line.		✓
speech over the phone line.		Flawless
amily & Friends Mobile App Features	Flawless Basic	Plus
Maximum number of Family & Friends who can use the App	3	Unlimited
ocation map with precise location via Wifi, GPS and LBS	✓	✓
ocation history, date & time	✓	✓
nstant location refresh (request current location from Wristband)	✓	✓
/ristband status – Battery level	✓	✓
/ristband status – Current signal level	✓	✓
/ristband status – Currently connected/not connected	✓	✓
/ristband status – On charge/being worn/not worn		✓
Vristband status – Connection quality – indication of mobile signal strength last 24 hours		✓
√ristband status – Activity level & sedentary alerts		✓
Vristband status – Advanced battery use/charging graph to identify low battery & battery mpty periods.		✓
vent history – summary of events	✓	✓

Alerts history – summary of important alerts such as fall detected or critical battery	✓	✓
Family & Friends chat messaging between App users.	✓	✓
Health measurements – Heart rate (beats per minute)**		✓
Health measurements – Blood pressure (systolic mmHg & diastolic mmHg)**		✓
Health measurements – Body temperature (degrees Celsius)**		✓
Health measurements – Blood oxygen %**		✓

^{**} Important Note About Health Measurements – The measurements given by your Wristband are Indicative Only and should not be relied upon for professional or medical purposes.

	1

	1